


CATALOGUE

English





History of DENTECH

In 1907, Dentech Corporation, a pioneering manufacturer in Dental Instruments, was founded in Yanaka, Tokyo by Kenjiro Kimura, and it was originally named as Kimura Kousakusho.

We manufactured the first Japan made dental burs.

After Second World War

In 1944, the first factory was destroyed due to Second World War. In 1947, after the war, we transferred our factory to Itabashi in Tokyo and reconstructed it.

However, in 1954, the strongest ever typhoon unfortunately struck around the area and the factory collapsed again.

In 1955, we rebuilt a new factory at the place nearby the destroyed factory to continue providing our products for the customers.

First Advancement Overseas

After Second World War, a lot of dentists came to Japan from the U.S. accompanying the occupation army and needed trustworthy dental instruments with high quality.

Thus we were able to develop our products through the cooperation of them.

In 1971, we launched to export our products overseas and our first overseas customer was in the United States.

Renamed as DENTECH

Kimura Kosakusho was renamed as Dentech Corporation in 1974. Since 1976 we have manufactured drills for implantology as the first implant drill manufacturer in Japan.

Nowadays most Japanese implant makers adopt our drills.

The Current DENTECH

Today Dentech is well known as one of the companies which have the longest over 100 - year history in dental industry.

KSK rubber dam instruments, implantation drills and burs as well as many unique products with Japan-made high quality are well supported by dentists and proffesionals all over the world.



Joint Performance with Japanese

CRAFTSMANSHIP

and Latest Machine Technologies



We possess CNC lathe, CNC Grinding machine, Heat treating machine, 3D CAD/CAM systems, and other various latest machines and equipments at our own factory in Japan. At the same time we are proud of the greatest technique by our craftsmen whose skill is unequaled with any others. We have an advantage in manufacturing extremely delicate dental/medical instruments using these two wheels which are identified by the machine automation and craftsmen's tequnique.

We have been developing, designing and manufacturing any kinds of dental and surgical instruments needed in clinical site around the world with outstanding quality for over a century.

It is a great honor for us that our products are highly respected all over the world. Our products would surely satisfy you.

Dentech has been assessed and certified as meeting the requirements of ISO 13485.



KSK Rubberdam system's products has qualified for the European CE mark medical device directive.



CONTENTS

ENDODONTICS		P.6
KSK Rubber Dam Clamps	- P.6	Peeso Reamers P.22
Dental Dams	P.17	RTP Reamers P.22
Rubber Dam Accessories	P.18	Finger Ruler P.23
Rubber Dam Set	P.21	X-Ray Mesh Gauges P.23
M Peeso Reamers	P.22	Resin Post Methacrylics P.23
		Precision Root Canal Model PRO P.24
PROSTHETICS		P.26
Matrix Bands	P.26	Tooth Separators P.27
Matrix Retainer	P.26	Articulating Paper Holders P.28
Tooth Contact Gauges	P.27	New Fine Lock Trays P.29
New Inlay & Crown Setters		
IMPLANTOLOG	Υ	P.30
Twist Drills Black	P.30	Drill Extenders P.35
Step Drills For Guide	P.31	Drill Extension SHORT P.35
Twist Drills For Guide	P.31	Para Drill Aid II P.36
Twist Drills	P.32	Technika Bone Mill Set P.36
Trephine Burs	P.33	Technika Repair Tap P.37
Bone Trephines	P.33	Depth Gauges P.37
Circuler Knives	P.33	3.75 Special Drills P.37
Surgical Burs	P.34	Gingiva Gauge P.38
Implant Drill checker	P.35	Color Boxes P.38
		Bur Stands P.38
LABORATORY -		P.39
Vulcanite Burs	P.39	Air Jet Diamonds P.39
MISCELLANEO	US -	P.40
	P.40	Parallel ReflectorP.42
Dry Field Mouth Props	P.40	pinlight2P.43
Lip Expanders	P.40	Excavators ————————————————————————————————————
Bite Expander	P.41	Wax Spatulas 7A P.44
Big Oral Mirrors Real-Fine	P.41	Evans P.45
Cotton Pliers	P.42	Hollenback #3 P.45
Tongue Guards & Saliva Ejector	P.42	Woodson #1 P.45

-Autoclavable-










KSK Rubber Dam Clamps






Features


- KSK Rubber Dam Clamps are made of specially selected and tempered stainless steel.
- Ideal spring action and high resistance against the corrosion.
- Various forms are available, so you can chosse clamps in accordance with the teeth.
- 4000 times open-close metal fatigue test goes for every production lot to assure long-term spring action.
- For stable forms, our craftsmen adjust bending every KSK clamps one by one manually.
- All beaks are sharpened for best fitting on the teeth.
- "BLACK" and "MATT" can avoid eye fatigue of dentists from reflecting operation light.




[Important Safety Instructions]

It is likely to transform and to break when force is excessively added. Please do not expand rubber dam clamps more than the necessity.









beaks. Similar to 3 but wingless.

01012007

01012008




beaks.





For upper molars with festooned beaks.

For big lower molars with flat beaks.





R L

Shifted to 7A with similar form.


For big upper molars with festooned beaks. Same size as 7.

Molar Clamps


Product Code Normal Matt
Black

01012010

01014010









10

For upper molars with deeply festooned beaks which can grab subgingval part.

Similar to 8A but wingless.

For small lower right molars with flat beaks, paired with 11.






01012912

01012914

01014914





For small lower left molars with flat beaks, paired with 10.

Especially for the third upper left and lower right molars with flat serrated beaks, paired with 13A.

Especially for the third lower left and upper right molars with flat serrated beaks, paired with 12A.

Discontinued





R I

For partially erupted molars with deeply festooned beaks.

For partially erupted molars with deeply festooned beaks, bigger than 14.

Discontinued Shifted to 14A with similar form

Molar Clamps

01012015

Product Code




01012024


01012028

01012042

01014042

01014024






01012018

01012026

01014026





For big lower molars with flat and longer beaks, a strong spring.

For upper & lower molars with flat beak.

For upper left and lower right buccal molars' subgingival cavities. The buccal side has wider fringe slightly festooned beaks. Paired with 25.









For lower left and upper right buccal molars' subgingival cavities. The buccal side has wider fringe slightly festooned beaks. Paired with 24.

For upper molars in case of exposure of the tooth neck with slightly festooned beaks. Available for the second bicuspids.

For upper and lower molars with flat beaks, wide view.







For lower left molar subgingival with slightly festooned beaks. The projection on the buccal beak slips into the furcation. Paired with 31.


For lower right molar subgingival with slightly festooned beaks. The projection on the buccal beak slips into the furcation. Paired with 30.


For upper left molars with flat beaks, paired with 43.


Molar Clamps

Product Code Normal Matt Black

01012056








For upper right molars with flat beaks, paired with 42.

For buccolabial molar furcations with slightly festooned beaks. Enables to press down the gingiva.

01012067

For lower molars with with festooned beaks.






138 01012138 01014138 01019138


For big lower molars with very slightly festooned beaks. One size larger than 56.

Best for upper second molars with festooned beaks, compact size.

Discontinued

For upper right and lower left third molars with serrated offset beaks, paired with 139.





201 01012201 01014201 01019201

For upper left and lower right third molars with serrated offset beaks, paired with 138.

Discontinued Shifted to 7 with same form.

For upper molars with slightly festooned beaks.


Molar Clamps

Product Code




01012204

50580003




For big lower molars with flat beaks. A little bigger than 200.







For lower right molars and the first bicuspids with flat beaks, Paired with 204.




For lower left molars and the first bicuspids with flat beaks, Paired with 203.




For upper molars with very slightly festooned beaks.





For lower left molars with cheek and tougue pressure shoving wings. No need to wear dental dam but gauze or cotton rolls. Paired with G-2.



For lower right molars with cheek and tougue pressure shoving wings. No need to wear dental dam but gauze or cotton rolls. Paired with G-1.







For small bicuspid with flat beaks. Available for deciduous anteriors as well.

Bicuspid Clamps

Product Code Normal Matt
Black

01012002


01014002


01019002

01012022

01014022

01012206





01012901

01012992

01014992


01012029


R L

For bicuspid root with festooned beaks.

For cervical cavities of bicuspid with festooned beaks.

For lower bicuspid with flat beaks.






01019022

For bigger upper and lower bicuspid with flat beaks. A little bigger than 2.

Similar to 2A but wingless.

For upper and lower bicuspid with flat beaks, similar to 207 but wingless.







For upper and lower bicuspid with festooned beaks, similar to 206 but wingless.


For upper and lower bicuspid with flat beaks, for wider view.

For upper and lower bicuspid with festooned beaks. Similar to 27 but winged.

Bicuspid Clamps

Product Code






208








01012208

01014208


01019208





For upper and lower bicuspids with flat beaks. Similar to 22 but winged.

For big bicuspid. Similar to 207 but the fringe is little bigger.

For lower bicuspid with flat beaks. Similar to 208 but a little smaller and stronger spring.










For anteriors with slightly offset beaks.

For labial cavities on anteriors with slightly offset beaks.








Offset beaks and stronger spring.


For anteriors labial and canines and cavities on bicuspid labial .


Labial Clamps

Discontinued

Normal **Product Code** Matt Black







Discontinued Shifted to 9 with same form.

For gum compression on anteriors with offset beaks.


Similar to 212 but a little stronger pressure to gingiva.

* #214 is for sale while stock lasts.












50050000

For exposure of the tooth neck and for the temporary retraction of the gingiva. Screw fixation guarantees a film grip. Ideal e.g. for adhesive fixing of veneers.







01013005

Deciduous Clamps

Product Code





















R L



For upper left the first deciduous molars. (I/64)

For upper left the second deciduous molars. (J/65)


For lower left the first deciduous molars. (L/74)

















For lower left the second deciduous molars. (K/75)

For lower right the first deciduous molars. (S/84)


For lower right the second deciduous molars. (T/85)

















For lower left the first deciduous molars. (L/74) Paired with P-2.

For lower right the first deciduous molars. (S/84)
Paired with P-1.








serrated beaks.

For bicuspids, especially for subgingival cavities with deeply festooned serrated beaks.





For lower bicuspids with flat ser-


rated beaks.












Used for labial cavities on anteriors with offset serrated beaks.





For lower molars with festooned beaks.

Dental Dams




Dental Dams (Pure natural rubber)

Features

- KSK dental dams have a good resistance and elastic property.
- Packed in airtight aluminium bags to keep the rubber sheet fresh.
- "Thin" is made of color free material then capable for monitoring the condition of an patient during isolation treatment.







Dental Dam Non-Latex

Features

- Used when there is a concern regarding an individual's possible allergic reaction to the proteins present in natural rubber.
- Packed in airtight aluminium bags to keep the rubber sheet fresh.

[Specification]

Packaging: 20 sheets / pack | Size: 6" x 6" (152mm x 152mm)

Thickness: 0.20mm | Color: Blue | Scent: Non

Product Code 01023040
Product Dental Dam Medium (cut)6"x6" Non-Latex



Dental Dam Non-Latex Powder-Free

Features

- Powder-free minimizes the risk of illness or injury caused by allergy for all dental health care personnel and patients.
- Packed in airtight aluminium bags to keep the rubber sheet fresh.

[Specification]


 $\begin{aligned} & Packaging: 15 \ sheets \ / \ pack & Size: 6" \ x \ 6" (152mm \ x \ 152mm) \end{aligned} \\ & Thickness: 0.25mm & Color: Purple & Scent: Mint \end{aligned}$

Product Code 01023050

Product Dental Dam Heavy (cut)6"x6" Non-Latex Powder-Free



Rubber Dam Accessories









eature:

Original design of DENTECH Rubber Dam Clamps are automatically off from the forceps at the place of operator fits them on tooth.

[Specification]

Material: Stainless steel
Size: 165 x 70mm

Product Code 01031040

Product KSK Rubber Dam Clamp Forceps







Rubber Dam Frames Young Type









Rubber Dam Template







Rubber Dam Clamp Boards









[Specification]

Material: Stainless steel Size: 90 x 90mm

With 8 clamps: #9 x 2, 2A, W2A, 5, 56, 203 & 204 or Board only (without clamps)

Product Code 01041000
Product Rubber Dam Clamp Board No.8 without Clamps
Product Code 01040000
Product Rubber Dam Clamp Board No.8 with Clamps
Product Code 01040080
Product Code 01040080 Product Rubber Dam Clamp Board No.8 with Clamps "Matt"



No.12

*From October 1, 2018, #211 was shifted to the same form #9.

[Specification]-

Material: Stainless steel Size: 150 x 90mm

With 12 clamps: #210, 9, 0, 2, 2A, 206, 207, 5, 56, 7, 201 & 202 or

Board only (without clamps)

board only (menode clamps)	
Product Code 01041010	
Product Rubber Dam Clamp Board No.12 without Clamp	s
Product Code 01040010	
Product Rubber Dam Clamp Board No.12 with Clamps	
Product Code 01040090	
Product Rubber Dam Clamp Board No.12 with Clamps "Matt	."
Product Code 01040060	
Product Rubber Dam Clamp Board No.12 with clamps "Black	"



No.18

*From October 1, 2018, #211 & #200 were shifted to the same form #9 & #7 respectively.

[Specification]

 $Material: Stainless \, steel \quad \, Size: 160\,x\,120mm$

With 18 clamps : #210, 9, P-1, P-2, 22, 27, 29, 208, 209, 18, 138, 139, 7,201,202, 203, 204 & 205 or Board only (without clamps)



Rubber Dam Clamp Boards







Deciduous Clamp Board

Material: Stainless steel Size: 90 x 90mm Complete kit of 8 patterns with special Clamp board. #55-upper right E #54-upper right D #64-upper left D #65-upper left E #85-lower right E #84-lower right D #74-lower left D #75-lower left E

v v
Product Code 01041061
Product Rubber Dam Clamp Board without Deciduos Clamps
Product Code 01041060
Product Rubber Dam Clamp Board with Deciduos Clamps
Product Code 01041070
Product Rubber Dam Clamp Board with Deciduos "Matt"
Product Code 01041080
Product Rubber Dam Clamp Board with Deciduos "Black"



I-Series

* I-Series is for sale while stock lasts.

[Specification]

Material: Stainless steel Size: 260 x 165mm

W8A, 10, 11, 12A, 13A, 14, 14A, 56, 56A (Total 28 Pcs.)

Product Code 01043010 Product Rubber Dam Clamp Board I-series without Clamps
Product Code 01043001 Product Rubber Dam Clamp Board I-series
Product Code 01043005 Product Rubber Dam Clamp Board I-series "Matt"
Product Code 01043003 Product Rubber Dam Clamp Board I-series "Black"



S-Series \star S-Series is for sale while stock lasts.

[Specification]-

Material: Stainless steel Size: 260 x 165mm

#210, 211, 212, 15, 18, 22, 24, 25, 26, 27, 28, 29, 30, 31, 51, 138, 139, 200, 201, 202, 203, 204, 205, 206, 207, 208, 209, G-l, G-2, P-l, P-2 (Total 31 Pcs.)

Product Code Product	01043020 Rubber Dam Clamp Board S-series without Clamps
Product Code Product	01043002 Rubber Dam Clamp Board S-series
Product Code Product	01043006 Rubber Dam Clamp Board S-series "Matt"
Product Code Product	01043004 Rubber Dam Clamp Board S-series "Black"



Rubber Dam Set


Rubber Dam Technic Set

Contents: 1) Clamp Bord #12

- ②Rubber Dam Frame Large(Adult)
- ③Rubber Dam Punch(Ainsworth)
- $@{\tt Rubber\,Dam\,Clamp\,Forceps}$

Clamps: #210, #9, #0, #2, #2A, #206, #207, #5, #56, #7, #201, #202We have 3 colors of rubber dam clamps. Please inquire at the time of your order.





*From October 1,2018, #211 was shifted to #9 with same form.





Rubber Dam Starter Kit

[Specification]-


Basic assorment for rubberdam technique

- Contents: 1 Clamp Board #12 with clamps
 - ②Rubber Dam Frame(Adult)
 - ③KSK Rubber Dam Punch(Ainsworth)
 - **4** KSK Rubber Dam Clamp Forceps
 - ⑤ Dental Dam Non-Latex Blue(20 sheets)
 - **©** Rubber Dam Template

Product Code 01030130

Rubber Dam Starter Kit





*From October 1,2018, #211 was shifted to #9 with same form.






M Peeso Reamers

Material: Stainless steel
Sharp shaped top
Pilot hole drills for root canal treatment.

Product Code 02312200 Produc	M Peeso Reamer CA #000 (set of 4 pcs)
Product Code 02312201 Produc	M Peeso Reamer CA #00 (set of 4 pcs)
Product Code 02312202 Produc	M Peeso Reamer CA #0 (set of 4 pcs)
Product Code 02312203 Produc	M Peeso Reamer CA #1 (set of 4 pcs)
Product Code 02312204 Produc	M Peeso Reamer CA #2 (set of 4 pcs)
Product Code 02312205 Produc	M Peeso Reamer CA #3 (set of 4 pcs)




Peeso Reamers

[Specification]

Material: Stainless steel

Round shaped top

Pilot hole drills for root canal treatment.



RTP Reamers

[Specification]

Material: Stainless steel

Round shaped

Final drills for root canal treatment.





			mm
	L	φD	arphid
# 1	9	1.35	0.92
# 2	10	1.65	1.17
# 3	10	1.95	1.48

#3

Finger Ruler

Can be attached on finger by free size. The adjustable stopper keeps it's position by spring.

[Specification]-

Material: Stainless steel Measurement length: 32mm

Product Code 03050000 Product

Finger Ruler







X-Ray Mesh Gauge

Description

Suitable for root canal treatment. Accurate examination before placing an implant. checking of calculus quantity and its area.

[Specification]-

Cover Material: Polyestel Mesh Line Material: Stainless steeel

- Dental Size
 - Size: 40 x 30mm Fine Grid: pitch: mm Wide Grid: pitch 5mm

Size:300 x 150mm Fine Grid: pitch 2mm Wide Grid: pitch 10mm


Product Code 07050010 Product

X-Ray Mesh Gauge, Dental

Product Code 07050011

Product X-Ray Mesh Gauge, Panolama





Resin Post Methacrylics

Description


You use the core build upon by the indirect method. It disappears, becomes gas when casting. For the post of the same size and RTP reamer, post selection is easy.

[Specification]-

Fine plastic burn - out posts.

Indirect only

Product Code	05000050
Product	Resin Post Methacrylic #R1 (set of 50 pcs)
Product Code	05000051
Product	Resin Post Methacrylic #R2 (set of 50 pcs)
Product Code	05000052
Product	Resin Post Methacrylic #R3 (set of 50 pcs)



Precision Root Canal Model PRO MTO

Features

- Reproduced precisely base on human root canal.
- For practical training;

Line Up: 24 types, 3 Levels (Basic • Normal • High) by difficulty level

Having visible root shape.

[Specification]-

 $Main\ Body: Transparent\ resin$ Root Part: Wax (white) filling

MOQ: 150 pcs./box (any type of combination is available)

Product Code As mentioned below As mentioned below



Inventor: Dr. Takatomo Yoshioka


- · Graduated from Graduate school of Tokyo Medical and Dental Univ. in 1996.
- By his lots of experiences of Conservative therapy and Endodontic Therapy awarded
- Excellent Treatise award, Gramd, $Grand\ chairman's\ award,\ Academic\ award.$
- After opening Endodontic trearment clinic at Tokyo in 2010, also focusing on training new dentists with his case study meeting.

Level		Product Code	Product		Classification						
	•	30000000	LIW2-01-01	Mandibular anterior teeth	Туре	Π	Orifice	2	Apical foramen	1	
	8	30000010	LM6W3-03-01	Mandibular molar teeth	Туре	III	Orifice	2	Apical foramen	2	
BASIC	8	30000020	LPMW1-03	Mandibular premolar teeth	Туре	I	Orifice	Γ	Apical foramen	1	
DASIC	4	30000030	U1-04	Maxillary anterior teeth	Туре	I					with lateral branch
	6	30000040	UMW3-03-2	Maxillary molar	Туре	III	Orifice		Apical foramen	2	medullary space opened version of UMW3-03
	6	30000050	UPMW2-07	Maxillary premolar	Туре	П	Orifice		Apical foramen	1	
	7	30000060	UPMW3-05-3	Maxillary premolar	Туре	III	Orifice	2	Apical foramen	2	
	8	30000070	LIW1-03	Mandibular anterior teeth	Туре	Ι	Orifice		Apical foramen	1	
	9	30000080	LM6W3-01	Mandibular molar teeth	Туре	III	Orifice	2	Apical foramen	2	"distal root apical foramen curved rapidly to distal"
	1	30000090	LMW1-02	Mandibular molar teeth	Туре	I	Orifice		Apical foramen	1	apical foramen curved gently
	①	30000100	LPMW1-06	Mandibular premolar teeth	Туре	I	Orifice		Apical foramen	1	
NORMAL	®	30000110	U1-05	Maxillary anterior teeth				Γ			
	B	30000120	UMW3-01-2	Maxillary first molar	Туре	III	Orifice	Γ	Apical foramen	2	narrower root canal than UMW3-01
	1	30000130	UPMW2-12	Maxillary premolar	Туре	II	Orifice		Apical foramen	1	
	(b)	30000140	UPMW3-01	Maxillary premolar	Туре	III	Orifice	2	Apical foramen	2	
	16	30000150	LIW1-01	Mandibular anterior teeth	Туре	I	Orifice	T	Apical foramen	1	
	1	30000160	LM6W2-04	Mandibular molar teeth				Г		Γ	
	®	30000170	LM6W3-02	Mandibular molar teeth	Туре	Ш	Orifice	2	Apical foramen	2	with distolingual root
	1	30000180	LPMW4-01	Mandibular premolar teeth	Туре	IV	Orifice		Apical foramen	2	divided buccolingually in the middle of the root
HIGH	20	30000190	U1-01	Maxillary anterior teeth	Туре	I					central incisor, without lateral branch
HIGH	a	30000200	UMW3-02	Maxillary second molar	Туре	III	Orifice	Γ	Apical foramen	2	coalesced MB root and P root including thin root,
											seen in upper second molars
	22	30000210	UPM-3roots-01	Maxillary premolar	Туре	V	oth	ner :	forms		3 roots including 2 roots in buccal **
	23	30000220	UPMW2-24-C	Maxillary premolar	Туре	П	Orifice		Apical foramen	1	with caries **with lateral branch
	24	30000230	UPMW4-02	Maxillary premolar	Туре	IV	Orifice		Apical foramen	2	

BASIC :For Dental students and interns

NORMAL :For Graduated New Dentists

HIGH: For Experts



PROSTHETICS

Matrix Bands

132°

Matrix Bands

[Specification]		
• •	ness : 0.04mm	
Product Code 04021101	Product	Matrix Bands #A (set of 100 pcs)
Product Code 04021102	Product	Matrix Bands #B (set of 100 pcs)
Product Code 04021103	Product	Matrix Bands #C (set of 10 pcs)
Product Code 04021104	Product	Matrix Bands #D (set of 10 pcs)
Product Code 04021105	Product	Matrix Bands #E (set of 10 pcs)
Product Code 04021106	Product	Matrix Bands #F (set of 10 pcs)
Product Code 04021107	Product	Matrix Bands #G (set of 10 pcs)
Product Code 04021108	Product	Matrix Bands #H (set of 10 pcs)
Product Code 04021109	Product	Matrix Bands #I (set of 100 pcs)
Product Code 04021110	Product	Matrix Bands #J (set of 100 pcs)


#A:General		#F	
#B		#G	3500
#C		#H	
#D	#I:Ge	enera • C	hild
#E		#J:Child	
	TYPE	Length	Width
	PARTITION	Length 63	Width
	TYPE A B	Length 63 63	-
	A	63	Width 6
Length	A	63 63	6 7.5 8
Length Width	8 C	63 63 65	6 7.6
k	B C D	63 63 65 51	Width 6 7.6 8 5
k	A B C D	63 63 65 51 66	Width 6 7.6 8 5
K	A B C D	63 63 65 51 66 64	Width 6 7.6 8 5 5
K	A B C D F G	63 63 65 51 66 64 49	Width 6 7.6 8 5 5 5 4.6


T Bands

Description

Can be used without the Matrix Retainer.







Matrix Retainer #9-T Tofflemire



Description

By a screw that has been improved, it will be a fine movement smooth.

[Specification]

Material: Stainless Steel Size: 64mm

Product Code	04019200
Product	Matrix Retainers Tofflemier Type, Universal Adult, 9-T



Tooth Contact Gauges

132°C

Silicone rubber : Dry heat are prohibited

Description

For checking interdental contact points. In three colour coded thicknesses.

[Specification]-

Material: Stainless Steel Holder size:130 x 22 mm

Thickness: Blue 0.05 mm, Yellow 0.11 mm, Red 0.15 mm

Set: Holder x lpc, Blue Gauges x 3pcs, Yellow Gauge x lpc, Red Gauge x lpc





Product Code	50050022
Product	Tooth Contact, Refill Gauge Yellow (0.11mm)
Product Code	50050023

Tooth Contact, Refill Gauge Red (0.15mm)

New Inlay & Crown Setters

[Specification]-

 $Material: Stainless\ steel\ and\ silicone\ rubber$

Long handle for molars has short pushing pin(SIZE : 79 x 24mm)

Short handle for bicuspids has long pushing pin (SIZE : $69 \times 24 \text{mm}$)

Product Code 06092300
Product New Inlay & Crown Setter Set
Product Code 06092310
Product New Inlay & Crown Setter L Set
Product Code 06092320
Product New Inlay & Crown Setter S Set
Product Code 50050035
Product New Inlay & Crown Setter, Extra Rubber Pads 2 pcs



Tooth Separators

Elliot's Separator

Description

The pushing screw can be attached both right and left side. Improved jaw design enable to be stable while pushing by the screw.

[Specification]-

 $Material: Stainless\ Steel$

Size: 37 x 67mm

Product Code	01071000
Product	Elliot's Separator



PROSTHETICS

Tooth Separators

Ivory's Separator

With improved curved jaw design, thus separators become stable while pushing by the screw.

It is applied in the front area of the upper and lower teeth.

[Specification]-

Material: Stainless Steel Size: 37 x 67mm

roduct Code 01072000

Ivory's Separator



Articulating Paper Holders

Articulating Paper Holder, M-Type

It is made with milled in order not to move articulating paper and designed to be convenient for the operator.

[Specification]

Material: Stainless Steel

Size: 145mm

roduct Code **06121000**

Articulating Paper Holder, M-Type



Articulating Paper Holder, Y-Type

Sure holding of Articulating papers by pins and holes. Ideal to make for full dentures.

[Specification]-

Material: Stainless Steel

Size: Large 74 x 60mm Small 66 x 48mm

Product Code Product	06122000 Articulating Paper Holder, Y-Type, Large
Product Code Product	06122100 Articulating Paper Holder, Y-Type, Small



New Fine Lock Trays

Description

"Retiention rims" types. Brass , chrome plated . Complete with set of $8\,\mathrm{size}$ or individually each size .

Upper				mm
	U11	U12	U13	U14
Wide	8 5	8 1	7 5	6 8
Depth	5 9	5 7	5 6	5 4
Height	19	2 0	2 0	2 1

Lower	mm			
	L11	L12	L13	L14
Wide	8 0	7 6	7 6	7 2
Depth	6 3	6 5	6 1	6 2
Height	2 1	2 1	2 2	2 3



Product Code 06040000
Product New Fine Lock Tray Assort (set of 8 pcs)
Product Code 06041011
Product New Fine Lock Tray #U11
Product Code 06041012
Product New Fine Lock Tray #U12
Product Code 06041013
Product New Fine Lock Tray #U13
Product Code 06041014
Product New Fine Lock Tray #U14
Product Code 06042011
Product New Fine Lock Tray #L11
Product Code 06042012
Product New Fine Lock Tray #L12
Product Code 06042013
Product New Fine Lock Tray #L13
Product Code 06042014
Product New Fine Lock Tray #L14



















IMPLANTOLOGY

Twist Drills Black

Description

A surgical drill with increased abrasion resistance, corrosion resistance, durability and hardness with DLC coating.

It comes with a movable stopper which allows for the length of the drill to be adjusted and a dedicated wrench. The adjustment range for the stopper is between 6 - 11mm for short, and 11 - 16mm for long.

What's DLC?

DLC is short for "Diamond Like Carbon", and it is technology that coats a nano - level hard film with properties that resemble that of diamonds onto metallic surfaces. While having approximately 5 times the hardness compared to existing products, it has shown overwhelming low friction factor and it is recognized as the most technologically advanced among the many surface treatment techniques available.


[Specification]-

Material: Stainless steel


Product Code	30073000	Product	Twist Drill "Black" Short 2.0S
Product Code	30073010	Product	Twist Drill "Black" Short 2.6S
Product Code	30073020	Product	Twist Drill "Black" Short 2.8S
Product Code	30073030	Product	Twist Drill "Black" Short 2.9S
Product Code	30073040	Product	Twist Drill "Black" Short 3.0S
Product Code	30073050	Product	Twist Drill "Black" Short 3.1S
Product Code	30073060	Product	Twist Drill "Black" Short 3.2S
Product Code	30073070	Product	Twist Drill "Black" Short 3.3S
Product Code	30073080	Product	Twist Drill "Black" Short 3.4S
Product Code	30073090	Product	Twist Drill "Black" Short 3.5S
Product Code	30073100	Product	Twist Drill "Black" Short 3.6S
Product Code	30073110	Product	Twist Drill "Black" Short 3.7S
Product Code	30073120	Product	Twist Drill "Black" Short 3.8S
Product Code	30073130	Product	Twist Drill "Black" Short 4.0S
Product Code	30073140	Product	Twist Drill "Black" Short 4.2S
Product Code	30073150	Product	Twist Drill "Black" Short 4.5S
Product Code	30073160	Product	Twist Drill "Black" Long 2.0L
Product Code	30073170	Product	Twist Drill "Black" Long 2.8L
Product Code	30073180	Product	Twist Drill "Black" Long 3.1L
Product Code	30073190	Product	Twist Drill "Black" Long 3.2L
Product Code	30073200	Product	Twist Drill "Black" Long 3.3L
Product Code	30073210	Product	Twist Drill "Black" Long 3.5L
Product Code	30073220	Product	Twist Drill "Black" Long 3.8L

Short

2.05 2.65 2.85 2.95 3.05 3.15 3.25 3.35



Long




Step Drills For Guide

132° ⟨⟩⟨

Description

Simplification of drill step has been made possible. It is ideal for pilot holes for tapered implants.



Twist Drills For Guide

Description

When using the 10mm surgical guide, it is marked to drill depths of 7mm, 8.5mm, 10mm.

[Specification]

Material: Stainless steel

Product Code 30074020	
Product Twist Drill "Black" for guide φ2.0mm	
Product Code 30074030	
Product Twist Drill "Black" for guide φ3.0mm	



IMPLANTOLOGY

Twist Drills



Description

Our drills are made of stainless steel with high carbon. It assures enough hardness, good cutting efficiency and long life. But rather weak for rust. Careful maintenance will be required to avoid rust and to assure along life.

Good maintenance

- $\textcircled{1} \ Immediately after use, drills should be dipped into anti-rust solution. }$
- ② Stopper should be removed while cleaning and autoclave.
- 3 Blood and bones should be removed completely after use.
- ④ Use only refined water for autoclave.
- (5) Use drill stand while autoclave, do not put into sterilization bag directly.

Product Code 30030200	Product Twist Drill Short 2.0S
Product Code 30030280	Product Twist Drill Short 2.8S
Product Code 30030310	Product Twist Drill Short 3.1S
Product Code 30030320	Product Twist Drill Short 3.2S
Product Code 30030330	Product Twist Drill Short 3.3S
Product Code 30030350	Product Twist Drill Short 3.5S
Product Code 30030380	Product Twist Drill Short 3.8S
Product Code 30030400	Product Twist Drill Short 4.0S
Product Code 30030420	Product Twist Drill Short 4.2S

	205 285 315 325 335 355 385 405 425 455 555.
0	
	20L 28L 31L 32L 33L 35L 38L

Product Code 30030450	Product Twist Drill S	Short 4.5S
Product Code 30030550	Product Twist Drill S	Short 5.5S
Product Code 30030201	Product Twist Drill	Long 2.0L
Product Code 30030281	Product Twist Drill	Long 2.8L
Product Code 30030311	Product Twist Drill	Long 3.1L
Product Code 30030321	Product Twist Drill	Long 3.2L
Product Code 30030331	Product Twist Drill	Long 3.3L
Product Code 30030351	Product Twist Drill	Long 3.5L
Product Code 30030381	Product Twist Drill	Long 3.8L

Twist Drills Internal Cooling


Description

The drills with internal cooling system should be used also with external cooling. Drill speed should be controlled from 600rpm(thick drills) up to 1000rpm(for fine drills).


Cutting efficiency is coming bad by each use and bad drill will cause heat damages against bone. Only good drills should be used for operation. For very hard bones, it is recommended to use surgical burs prior to implant drills. Surgical burs will help easy drilling and long life. The fix screw of drill stopper is so small and it can not be bear a strong force. Drills can be used also without stopper.

[Specification]—

Material : Stainless steel







Product Code 30031201	Product	Twist Drill Long 2.0 LI
Product Code 30031281	Product	Twist Drill Long 2.8 LI
Product Code 30031311	Product	Twist Drill Long 3.1 LI
Product Code 30031321	Product	Twist Drill Long 3.2 LI
Product Code 30031331	Product	Twist Drill Long 3.3 LI
Product Code 30031351	Product	Twist Drill Long 3.5 LI
Product Code 30031381	Product	Twist Drill Long 3.8 LI
Product Code 30032000	Product	Screw for Implant Drill Stopper
Product Code 30032001		
Product Hexagonal	Wrench for	Implant Drill Stopper

Trephine Burs

Trephine Burs do not have blade at center, thus it will be rather difficult to get smooth cutting for a start. Enough exercise is recommended before use. After Trephine Burs penetrate into bone, then smooth cutting is available. To avoid choking on bone inside Trephine Bur, pull of Trephine Burs up and down, so that bone will always remain at patient. Remained bone can be collected by surgical curettes. Enough cooling and possibly low speed is recommended to avoid heat damage against bone.

[Specification]-Material: Stainless steel

Product Code 30060240	Product	Trephine Bur 2.4
Product Code 30060270	Product	Trephine Bur 2.7
Product Code 30060300	Product	Trephine Bur 3.0
Product Code 30060330	Product	Trephine Bur 3.3



Product Code 30060380	Product Trephine Bur 3.8
Product Code 30060430	Product Trephine Bur 4.3
Product Code 30060500	Product Trephine Bur 5.0
Product Code 30060600	Product Trephine Bur 6.0

Bone Trephines

Bone Trephine is used for collection of big volume of bone. Enough exercise is recommended before use.

[Specification]

Material: Stainless steel

Product Code 30070500	Product	Bone Trephine 5.0
Product Code 30070600	Product	Bone Trephine 6.0
Product Code 30070800	Product	Bone Trephine 8.0
Product Code 30071000	Product	Bone Trephine 10.0
Product Code 30071200	Product	Bone Trephine 12.0




Circular Knives


For cutting soft tissue. Use with very low speed or by finger. Complete with finger grip for each piece.

[Specification]-

Material: Stainless steel

Product Code 50010000
Product Circular Knife 3.3 without finger driver
Product Code 50010010
Product Circular Knife 3.5 without finger driver
Product Code 50010020
Product Circular Knife 4.0 without finger driver
Product Code 50010030
Product Circular Knife 5.0 without finger driver
Product Code 30051001
Product Finger Driver for Circular Knife





IMPLANTOLOGY

Surgical Burs

Surgical Bur External Cooling System

Description


Round Burs

Round burs for bone cutting, ruff cutting edges to prevent clogging. Full group blades ensures excellent cutting efficiency.

Dr. Lindemann Bur

Dr. Lindemann Bur has fully around blades.

Product Code 30011400	Product Surgical Bur Round 141RF/010
Product Code 30011401	Product Surgical Bur Round 141RF/014
Product Code 30011402	Product Surgical Bur Round 141RF/018
Product Code 30011411	Product Surgical Bur Round 141RF/025
Product Code 30011412	Product Surgical Bur Round 141RF/031
Product Code 30011414	Product Surgical Bur Round 141RF/040
Product Code 30011415	Product Surgical Bur Round 141RF/050
Product Code 30011631	Product Dr. Lindemann Bur 163RF/014
Product Code 30011641	Product Dr. Lindemann Bur 164RF/018
Product Code 30011661	Product Dr. Lindemann Bur 166RF/021



Dr. Lindemann Bur




Surgical Bur Internal Cooling System

These hully group surgical burs with internal cooling system enable very smooth cut and prevent heat damages for bone. Cooling fluid flow out blood thus operator can get clear view of operating spot. Prepare #40 K Reamers to clean internal cooling hold.

Product Code 30014711
Product Surgical Bur Internal Cooling 471RF/023
Product Code 30011671
Product Surgical Bur Internal Cooling 167RF/023
Product Code 30014731
Product Surgical Bur Internal Cooling 473RF/031
Product Code 30014732
Product Surgical Bur Internal Cooling 473RF/040
Product Code 30014733
Product Surgical Bur Internal Cooling 473RF/050





Implant Drill Checker




Can be measured diameter of 2.0 - 5.9mm. Can be measured length of 6.0 - 18.0mm. With length scale of 5 - 20mm. Can be used for drill stand.

[Specification]-Material: Stainless steel

Product Code 30080012

Implant Drill Checker




Drill Extenders

Two types: External cooling/Internal cooling Extension length: External cooling is 16.5mm, Interna cooling is 18 mm.

[Specification]-

Material: Stainless steel







Product Code 30080001

Product Drill Extension External Cooling

Product Code 30080002 Product Drill Extension Internal Cooling

Drill Extension SHORT External Cooling



Can extend the drill length slightly 5mm. Easy attaching without tools.

[Specification]

Material: Stainless steel Diameter: 7.7 mm

Product Code 30079991

Drill Extension SHORT External Cooling



IMPLANTOLOGY

Para Drill Aid II

Description

- Attachment drill aid for parallel implantation.
- Four guid pins and three drill sleeves make surgery better and shorter.
- You can easily make a hole for parallel implantation.
- Pilot drill (dia 1.8 to 2.0mm) is available.
- All parts are autoclavable.

[Specification]-

Material: Stainless steel

Product Code 30080200

Para Drill Aid II





Technika Bone Mill Set

Description

- Simple and compact
- Easy operating

The angle mount enables easy and smooth operating.

Easy milling

The unique shaped cutter with DLC coating enable easy and fine milling.

Easy assembling

The simple structure saves time for operating and cleaning With sterization case and cleaning brush Cutter for substitute.

sell separately

Product Code Product	30080015 Technika Bone Mill
Product Code Product	30080016 Replacement Cutter
Product Code Product	30080117 Technika Bone Mill with sterilization case







Technika "REPAIR TAP"

- For internal threads repair.
- Made of high corrosion resistant stainless steel compatible with autoclave.
- Long lasting DLC coating.
- Safe Handle

[Specification]-

Material: Stainless steel (Tap and Handle) Surface Treatment : DLC coating (Tap) Thread Diameter x Pitch: M1.6 x 0.35 / M2.0 x 0.4

Made in Japan











Handle for Repair Tap

Compatible to both Tap

Depth Gauges

Graduations on 8, 10, 13, 15, 18, 20. Graduations can be exposed on X-ray film. Diameter: 1.3, 1.5, 1.9, 2.7, 3.0, 3.2, 3.4, 3.9mm.

[Specification]-

Material: Stainless steel






Product Code 30072300	Product Depth Gauge 3.0	
Product Code 30072320	Product Depth Gauge 3.2	
Product Code 30072340	Product Depth Gauge 3.4	
Product Code 30072390	Product Depth Gauge 3.9	

3.75 Special Drills

Diameter is 3.75mm. Twist drill with safety stopper.

Material: Stainless steel








Product Code	30043754
Product	φ3.75 Special Drill φ3.1 x 10
Product Code	30043755
Product	φ3.75 Special Drill φ3.1 x 13
Product Code	30043756
Product	$\phi 3.75$ Special Drill $\phi 3.1 \times 15$

IMPLANTOLOGY

Gingiva Gauge

The first division shows 1mm when top of measuring block and needle are coincide. Wide of graduation is 1mm and pitch is 2mm, thus each graduation shows the gingival thickness of 2, 3, 4, 5, 6, 7mm. Put the measuring block in right angle to gingiva, then prick gingival by light force until needle reaches to bone. Pull out carefully so that measuring block would not moves. Read the thickness at the back side of measuring block.

Product Code 30080013 Product Gingiva Gauge



- · Can be measured the thickness of gingiva from 2mm to 7mm in increments of 1mm during implant operation.
 •For a more exact and detailed measure. It is better to use
- micrometer caliper separately.

Color Boxes

- Sterilization Box Can be housed three bur stands.
- Color Boxes Box for small arrangement. The inside has been

[Specification]

Material: Aluminum

Size: 170 x 70 x 30mm (Sterilization Box) Size: 135 x 55 x 180mm(Color Boxes)

partitioned into six rows.

Color: Silver, Blue, Red, Gold, Green(Color Boxes)

Product Code 30080024	Product Sterilization Box
Product Code 30080018	Product Color Box "Red"
Product Code 30080019	Product Color Box "Green"
Product Code 30080020	Product Color Box "Blue"
Product Code 30080021	Product Color Box "Gold"
Product Code 30080022	Product Color Box "Silver"
Product Code 30080023	Product Color Boxes Set





Bur Stands

Description

You can stand 6 burs.

[Specification]-

Material: Stainless Steel

Bur Stand 40(w70 x D40 x H15mm) Compliant bur length: Max 35mm Bur Stand 45(W70 x D45 x HI5mm) Compliant bur length: Max 40mm

Product Code	30080010	Product	Bur Stand 40
Product Code	30080011	Product	Bur Stand 45





LABORATORY

Vulcanite Burs

[Specification]-

Material: Stainless Steel

Large and small types are available. Large: Maximum diameter is 6 mm. Small: Maximum dimeter is 4.8 mm.



* Straight hand piece type only.

Product Code 02101	.011 Product	Vulcanite Burs #A, LARGE
Product Code 02101	.012 Product	Vulcanite Burs #A, SMALL
Product Code 02101	.021 Product	Vulcanite Burs #B, LARGE
Product Code 02101	.022 Product	Vulcanite Burs #B, SMALL
Product Code 02101	.031 Product	Vulcanite Burs #C, LARGE
Product Code 02101	.032 Product	Vulcanite Burs #C, SMALL
Product Code 02101	.061 Product	Vulcanite Burs #F, LARGE
Product Code 02101	.062 Product	Vulcanite Burs #F, SMALL

Product Code 02111011	Product	Vulcanite Burs #AK, LARGE
Product Code 02111012	Product	Vulcanite Burs #AK, SMALL
Product Code 02111021	Product	Vulcanite Burs #BK, LARGE
Product Code 02111022	Product	Vulcanite Burs #BK, SMALL
Product Code 02111031	Product	Vulcanite Burs #CK, LARGE
Product Code 02111032	Product	Vulcanite Burs #CK, SMALL
Product Code 02111061	Product	Vulcanite Burs #FK, LARGE
Product Code 02111062	Product	Vulcanite Burs #FK, SMALL

Air Jet Diamonds

Features

Air Jet Diamond have air slits at cutting head, the slits supply strong air while cutting.

The strong air flow prevent the clogging of cutting materials.

Thus, you can very smoothly and effective cutting.

Product Code	02091001
Product	Air Jet Diamonds #1
Product Code	02091002
Product	Air Jet Diamonds #2
Product Code	02091003
Product	Air Jet Diamonds #3
Product Code	02091010
Product	Air Jet Diamonds #1A
Product Code	02091011
Product	Air Jet Diamonds #2A







MISCELLANEOUS

Lip & Cheek Retractor

Features

To retractor cheeks and to have a better view. Can be packed by quantities of customer's order.

[Specification]

Material : Polycarbonate

Material: Polycarbonate Set of 2pcs

Product Code 50570131 Product Lip & Cheek Retractor Set



Dry Field Mouth Props

Feature:

Great for air - abrasion prep and sealants.

- Comfortably holds the jaw open.
- Retracts and protects the patients cheek and tongue.
- Create a moisture free environment.

[Specification]

Material: Polycarbonate

Product Code 50570000	Product	Dryfield Mouth Props #1 labial (Set of 3 pcs)
Product Code 50570002	Product	Dryfield Mouth Props #2 large (Set of 3 pcs)
Product Code 50570001	Product	Dryfield Mouth Props #3 small (Set of 3 pcs)
Product Code 50570005	Product	Dryfield Mouth Props Assort (#1,2,3)



#1 Labial



#2 Large(Posterior)



#3 Small(Posterior)

Lip Expanders

Feature

Superb anterior access. Available in an adult, child or assorted (adult / child) set.

- Soft, plastic frame.
- Comfortable for the patient.
- Provides excellent anterior access.

[Specification]—

Material: Polypropylene

Product Code 50570010	Product	Lip Expander Adult Size (set of 2 pcs)
Product Code 50570011	Product	Lip Expander Child Size (set of 2 pcs)
Product Code 50570020	Product	Lip Expander Adult & Child, (set of 2 pcs)





Bite Expander

Description

"Bite Expander" has function of Retractor and Bite block. You can easily change 4 patterns for every treated area.

- This device serves the purpose of both a cheek retractor and bite block.
- Can be used on every part of the mouth.
- Made from silicone and polypropylene.

[Specification]

Material: Silicon rubber and Polypropylene

Product Code	07000150
Product	Bite Expander
Product Code	07000155
Product	Extra Rubber Pads for Bite Expander (2pcs)







Big Oral Mirrors Real-Fine

Features


- Ideal for distortion free intra oral photos.
- Extensive rhodium coating on one side.
- Sturdy, durable.
- With detachable silicone handle.
- Autoclavable in the provided fabric bag.

#1: Side views. #2: Occulusal views (medium). #3: Occulusal views (large). #4: For lingual side views of posteriors.

#5: For lingual side views of anteriors. Handle for #4, #5: aluminum.


Big Oral Mirror Set: Assort of #1,#2 and #3












MISCELLANEOUS



Tongue Guards & Saliva Ejector





escription

Three different sizes metal reflectors are changeable and they protect tongue from turbine.

[Specification]
Pipe: brass, chrohm plated Others: Stainless steel Size: 138 x 80mm





MOQ: 100 pcs. Please contact us about the price every order.

Parallel Reflector

Features

Mirror with parallel lines. Autoclavable.

[Specification]—

Coak type Material : Stainless steel

Size : \dot{H} 18 x W 48 mm, Thread Part M2.5 (P=0.45 L-5)

Product Code 50050100
Product Parallel Reflector



MOQ : 100 pcs

pinlight2

Features

- The mirror and the light guide are detachable and autoclavable.
- The handle is IP X7 waterproof and capable of alcohol disinfection.
- Much brighter 90lm White LED Lighting used for detecting caries.
- Roder's HR FRONT mirror comes with size 4.

[Specification]-

Light: 90LM, Cool white LED Mirror: HR FRONT coating, size 4

Battery : Two AAA-size

Waterproof Performance IPX7 JIS C 0920 (IEC 60529)

Continuous lighting time: 90 min

Product Code	30080350
Product	pinlight2
Product Code	30080420
Product	Mirror holder, set of 4pcs.
Product Code	30080360
Product	Light guide
Product Code	30080370
Product	Replacement O ring, set of 3pcs.



MOQ : 10 pcs

MISCELLANEOUS









Features

special hollow octagon handle which enables very light weight as well as ideal balance. high cutting performance and long life











MOQ : 100 pcs

■ Wax Spatulas #7A



Evans



Hollenback #3



Plastic Filling Instruments, Woodson #1

